

# See the difference with **Prosigna**<sup>®</sup>

## Breast Cancer Prognostic Gene Signature Assay



- PAM50-based report of intrinsic subtype
- The first genomic breast cancer assay that can be run in your local lab
- CE marked for FFPE tissue

## CE Product Summary Sheet

Prosigna is an in vitro diagnostic assay that uses the gene expression profile of cells found in breast tissue to assess a patient's risk of distant recurrence. The assay measures the gene expression profile using RNA extracted from Formalin-Fixed, Paraffin-Embedded (FFPE) breast tumour tissue. The gene expression data are weighted together with clinical variables

to generate both a subtype (luminal A, luminal B, Her2-enriched, or basal-like) and a score indicative of the probability of distant recurrence of disease. The assay is performed on the NanoString nCounter<sup>®</sup> Dx Analysis System using FFPE breast tumour tissue previously diagnosed as invasive breast carcinoma.

# Prosigna Sample Requirements




**Block  
Selected**



**H&E stain to identify  
tumour area and  
cellularity**



**Tumour area transposed  
to unstained slides and  
macrodissected**



**RNA extracted with  
manual kit**

The Prosigna assay is performed on RNA isolated from FFPE breast tumour tissue. A pathologist examines a hematoxylin and eosin (H&E) stained slide and identifies (and marks) the area of invasive breast carcinoma suitable for the test. The pathologist also measures the tumour surface area, which determines the number of unstained slides required for the test, and the tumour cellularity


to ensure the presence of sufficient tumour tissue for the test. A trained technician macrodissects the area on the unstained slides corresponding to the marked tumour area on the H&E-stained slide and isolates RNA from the tissue. The isolated RNA is then tested on the NanoString nCounter Dx Analysis System to provide test results including the Prosigna Score and risk category.

Specimen Attribute	Requirement	
Tissue input	Viable invasive breast carcinoma (ductal, lobular, mixed, or NOS)	
Tissue input format	Macro-dissected 10-micron-thick slide-mounted tissue sections	
Minimum tumour size	4 mm <sup>2</sup> tumour area	
Minimum tumour cellularity	10% within tumour area	
Minimum RNA amount	125 ng (12.5 ng/μL)	
Tumour area	≥ 100 mm <sup>2</sup> 20 – 99 mm <sup>2</sup> 4 – 19 mm <sup>2</sup>	1 slide 3 slides 6 slides

# Prosigna Algorithm and Output

The test is based on the PAM50 gene expression profile, and is performed on the nCounter Analysis system using RNA extracted from formalin fixed paraffin embedded (FFPE) breast tumour tissue samples. The algorithm uses a 50-gene expression profile to assign breast cancer to one of four PAM50 molecular subtypes determined by the tumour's molecular profile. The prototypical gene expression profiles (e.g. centroids) of the four PAM50 molecular subtypes were retrained on the nCounter Dx Analysis System using FFPE breast tumour samples collected from multiple clinical sites in North America. After performing the assay on a patient test sample, a computational algorithm based on a Pearson's correlation compares the normalized 50-gene expression profile of the patient test sample to the prototypical expression profiles of the four breast cancer intrinsic subtypes. The patient test sample is assigned the subtype with the highest Pearson's correlation.


The algorithm reports a Risk of Recurrence (ROR) score on a 0-100 scale, which is correlated with the probability of distant recurrence at ten years for post-menopausal women with hormone receptor-positive, early stage breast cancer. The report also provides a risk category (low, intermediate, or high). The ROR score is calculated using coefficients from a Cox model that includes the Pearson correlation to a 46-gene subset of the 50 genes<sup>1</sup> used to calculate each intrinsic subtype, a proliferation score, and tumour size. The test variables are multiplied by the corresponding coefficients from the Cox model to generate the score, which is then adjusted to a 0-100 scale based on coefficients generated from the training set of FFPE breast tumour samples. Risk categories are also reported based on cut-offs for ROR determined in a clinical validation study.



**Extract RNA from  
FFPE tumour sample**



**Run RNA and Prosigna CodeSet  
on nCounter Dx Analysis System**



**Patient specific  
expression profile**

<sup>1</sup>The 50-gene ROR and 46-gene ROR are functionally identical with a very high R-squared value equal to 0.997. Wallden et al. BMC Medical Genomics (2015) 8:54

# Ordering Information

## PRODUCT DESCRIPTION

## CATALOG NUMBER

## UNIT

### Prosigna® Gene Signature Assay

Complete kit for running Prosigna tests. Includes all CodeSet and Master Kit components; does not include RNA Isolation Kit.

PROSIGNA-001

PROSIGNA-002

PROSIGNA-003

PROSIGNA-004

PROSIGNA-010

One kit of 1 patient assay

One kit of 2 patient assays

One kit of 3 patient assays

One kit of 4 patient assays

One kit of 10 patient assays

### Roche FFPET RNA Isolation Kit

Includes 25 isolations per kit.

Roche-FFPET-025

Each

## CONTACT

email: [info@prosigna.com](mailto:info@prosigna.com)

For more information and details on how to offer Prosigna at your institution, please contact:

### EU Authorized Rep.

Emergo Europe  
Prinsessegracht 20  
2514 AP The Hague  
The Netherlands  
EU +31.70.345.8570

### NanoString Technologies, Inc.

530 Fairview Ave N  
Seattle, WA 98109  
Phone: 888-358-NANO  
(888-358-6266)



For more information, please visit [www.prosigna.com](http://www.prosigna.com)

Prosigna in conjunction with the nCounter® Dx Analysis System is CE marked for in vitro diagnostic use in post-menopausal women with Hormone Receptor-Positive (HR+), lymph node-negative, Stage I or II breast cancer and post-menopausal women with Hormone Receptor-Positive (HR+), lymph node positive (1-3 positive nodes, or 4 or more positive nodes), Stage II or IIIA breast cancer to be treated with adjuvant endocrine therapy. See Package Insert for further details at [www.prosigna.com](http://www.prosigna.com).

**nanoString**